Computer vision profiling of neurite outgrowth morphodynamic phenotypes

Ludovico Fusco^{1,7}, Fethallah Benmansour^{2,7}, Riwal Lefort^{3,7}, Kevin Smith^{4,7},
German Gonzales⁵, Catherina Barillari⁶, Bernd Rinn⁶, Pascal Fua², Francois Fleuret³
& Olivier Pertz¹

¹Institute for Biochemistry and Genetics, Dept. Biomedicine, University of Basel, Mattenstrasse 28, 4058 Basel, Switzerland

²Computer Vision Laboratory, Dept. I&C, École Polytechnique Féd'erale de Lausanne, Station 14, CH-1015 Lausanne, Switzerland

³Idiap Research Institute, Centre du Parc, Rue Marconi 19, 1920 Martigny, Switzerland

⁴Light Microscopy and Screening Center, Institute for Biochemistry, Schafmattstrasse 18, ETHZ-Honggerberg, CH-8093 Zurich, Switzerland

⁶Center for Information Sciences and Databases, Department of Biosystems Science and Engineering, Maulbeerstrasse, Building 1078, 1st Floor, 4058 Basel, Switzerland

⁷TODO

⁵TODO

Neurite outgrowth is a dynamic cell behavior that consists of morphogenetic processes such as neurite initiation, protrusion-retraction cycles, branching and growth cone navigation. Current knowledge of the underlying signaling events stem from molecular perturbations assessed at the steady state, which cannot capture the dynamic nature of this complex process. Here, we present NeuroDynamics an integrated pipeline to study neurite outgrowth dynamics. A microscopy platform allows for high content imaging of neurite outgrowth dynamics. A computer vision approach allows segmentation of cell shape and extraction of features describing neuronal morphodynamics. Statistical analysis then allows to automatically identify features that significantly discriminate between the control and perturbed state. We demonstrate the applicability of our approach by automatically annotating morphodynamic phenotypes in an RNA interference screen targeting a candidate Rho GT-Pase signaling network identified by a proteomics approach. NeuroDynamics is freely available as open source software to study neuronal outgrowth dynamics.

You are welcome to add some technical jargon given that we remain in the
 150 words limit!

• We can also think of another acronym than Neurodynamics!

Neurite outgrowth is essential to build the neuronal processes that connect the adult brain. This morphogenetic process is highly dynamic, and consists of a series of stochastic and repetitive events such as neurite initiation, elongation, branching, growth cone motility and collapse. Each of these different morphogenetic behaviors are likely to be regulated by different signaling networks. Furthermore, all these processes occur on length and time scales of microns and minutes to hours, suggesting an exquisite spatio-temporal regulation of the underlying regulating signals. This is not accessible with approaches in which the effects of molecular perturbations are assessed at the steady-state?,?, in which only snapshots of intrinsically dynamic behaviors are captured. Rather, understanding the signaling networks that regulate this complex process might benefit from understanding the effect of molecular perturbations on its dynamics. Advances in automated live cell microscopy now allows for high content image acquisition with high temporal resolution¹. However, from the analysis side, most computer vision approaches have until now only been used to automatically extract image features from cell populations at the steady-state 1,2, rather than from dynamic time-lapse datasets. In one study about mitosis, the temporal context in time-lapse datasets

has been taken into account to improve classification of functional cellular states by deconvolving ambiguous, transiently occurring phenotypes with similar morphology3. In another perturbation screen about cell scattering, global dynamic features such as cell migration speed have been measured4. However, no approaches have directly tried to automatically extract information about multiple dynamic cellular features to classify phenotypes. Here, we present Neurodynamics, an integrated pipeline that allows to study dynamic neurite outgrowth phenotypes in response to molecular perturbations. Neurodynamics consists of three entities: 1. a high content imaging platform for neurite outgrowth dynamics, 2. an automated neuronal shape segmentation algorithm that extracts morphological and morphodynamics features, 3. a feature selection algorithm that allows automatic inference of the neurite outgrowth morphodynamic phenotype. A scheme depicting the global approach is shown in Fig.1. We demonstrate the efficiency of our approach to study a variety of molecular perturbations.

Results

High content imaging of neurite outgrowth dynamics To study neurite outgrowth dynamics, we used the mouse N1E-115 cell neuroblastoma cell line, that can be easily triggered to extend neurites by simple serum starvation and plating

on the extracellular matrix protein laminin. To visualize soma and neurite morphology, we used a fusion of green fluorescent protein (GFP) with the F-actin binding peptide lifeact5. This provided high contrast on neurites and soma without excessive accumulation of fluorescence in the thick somata of differentiated N1E-115 cells, allowing for adequate imaging in N1E-115 cells by wide field microscopy (Figure 2a). For unambiguous cell identification, we also simultaneously labeled the nucleus with a mCherry fusion with a nucleus localization sequence (NLS), expressed from the same expression vector than Lifeact-GFP (Supplementary Fig.1a). We observed that transient expression of this construct did not affect neurite outgrowth in N1E-115 cells (Supplementary Fig.1b). To perturb different signaling pathways, we optimized a method to co-transfect the fluorescent marker and siRNAs simultaneously in cells, and observed efficient knockdown in response to knockdown of a panels of mRNAs (Supplementary note 1). To perform high content live cell imaging, we optimized our microscope for fast acquisition of multiple wells of a 24 well plate (Supplementary note 2). This allowed to acquire 240 fields of view in two channels across a 24-well plate with 12 minute time resolution for 20 hours. We performed different sets of experiments with 10 and 20x air objectives. Using 10 x objectives allowed for aquisition of a large field of view with typically 20 objects that could be consistently be observed for long time periods. The 20x objective allowed high resolution imaging of morphological features such as filopodia, but the moving cells often migrated out of the field of view resulting in loss of these objects. Typical 10x and 20x movies are shown (Supplementary movie 1 and 2). In these experiments, we tracked neuronal morphodynamics of differentiated cells that were replated on laminin, allowing to observe almost the whole neurite outgrowth process. We observed that illumination of neurons in the early phase of neurite outgrowth was toxic to the cells, leading to their death. We found that starting the observation process 3 hours postplating, a time at which cells also had extended 3-4 neurites, allowed the cells to survive until the end of the imaging process.

Describe verbally the neurite outgrowth process here?

Automated computer vision analysis of neurite ougrowth dynamics To capture the dynamic neurite outgrowth trajectories of neurons in the native and perturbed state, we developed a computer vision pipeline that allowed to automatically segment and track the soma and neurites in each frame of the timelapse datasets. Our approach first detects nuclei and associated somata at each time step. The nucleus of each neuron is detected as a Maximally Stable Extremal Region (MSER) ² from the Cherry channel. [As shown in the supplementary note

3, this is more robust than adaptive thresholding.] Using the detected nuclei as seed points, a region-growing algorithm segments the neurons soma. Next, the implemented multi-objects tracking algorithm ³ searches through the full set of nuclei and somata detections to extract the best K-shortest paths according to a similarity measure between two detections. The proposed similarity measure is derived from the Earth Movers Distance ⁴ between the intensity histograms of the detected somata regions (in the supplementary note 3, we show that this distance provides a good tradeoff between efficiency and precision). Neurons are detected and tracked at an accuracy of (TODO 95%) (still waiting for the GT to be completed for number crunching); Finally, the tracked somata are used to initialize a neurite segmentation and association algorithm based on shortest path computation and Voronoi tessellation. Comparison to manually annotated data demonstrates that (URGENT TODO: crunch the numbers: still waiting for the GT)

To Riwal

- 1. describe feature selection algorithms
- validation of feature selection algorithm by mixing data from two perturb states

- 3. explain approaches to take into account siRNA noise
- 1. Held, M. *et al.* CellCognition: time-resolved phenotype annotation in high-throughput live cell imaging. *Nature methods* **7**, 747–754 (2010).
- Nistér, D. & Stewénius, H. Linear time maximally stable extremal regions. In Proceedings of the 10th European Conference on Computer Vision: Part II, ECCV '08, 183–196 (Springer-Verlag, Berlin, Heidelberg, 2008).
- 3. Berclaz, J., Turetken, E., Fleuret, F. & Fua, P. Multiple object tracking using k-shortest paths optimization. *IEEE Transactions on Pattern Analysis and Machine Intelligence* (2011).
- 4. Pele, O. & Werman, M. Fast and robust earth mover's distances. In *ICCV* (2009).

Acknowledgements Put acknowledgements here.

Competing Interests The authors declare that they have no competing financial interests.

Correspondence Phone: +41 61 267 22 03; Fax: +41 61 267 35 66; E-mail: olivier.pertz@unibas.ch.